


ROOM PERIPHERALS

Patient handset



Description

Patient handset with JACK 6,5mm socket for connecting to LLC-50x call unit. It allows initiating a patient call, by pressing call push-button. Once a patient call has been started the green LED indicator stands steady ON until the call is cancelled by a nurse presence in the room. Easy to clean antimicrobial cover. Flush mounting.

SPECIFICATIONS

Socket (patient handset connector)	Angled male JACK 6,5mm
Power supply	24 V
Consumption (idle)	2 mA
Consumption (maximum)	6 mA
Dimensions (width x height x deep)	120x60x25 mm
Plastic	ABS / Polycarbonate
Colour (plastic)	Signal White (RAL 9003)
Operating Temperature	5°C to 55 °C
Protection index	IP 40
Mounting	Flush mounting (standard box)
Safety	EN 60065:2003
EMC	EN 50083-2:2007