



ROOM PERIPHERALS

Wall unit with call button



LLC-560



LLC-562

Description

Standard room peripheral with a patient call push-button with green LED indicator. It allows starting a patient call by pressing the call push-button or using connected patient handset. Once a patient call has been started the green LED indicator stands steady ON until the call is cancelled by a nurse presence in the room. Easy to clean antimicrobial cover. Flush mounting.

| SPECIFICATIONS | LLC-560 | LLC-562 |
|------------------------------------|------------------------------------------------------------------------|---------|
| Patient call button | Text: CALL | |
| | Icon: Nurse | |
| | Colour: Green (RAL 60118, yellow green) | |
| | Text / Icon colour: Black (RAL 9004, signal black) | |
| | Indicator led: Green | |
| Power supply | 24 V | |
| Consumption (idle) | 5 mA | |
| Consumption (maximum) | 20 mA | |
| Dimensions (width x height x deep) | LLC-560 (vertical): 75x120x15 mm LLC-562 (horizontal): 120x75x15 mm | |
| Plastic | ABS / Polycarbonate | |
| Colour (plastic) | White (RAL 9003, signal white) | |
| Operating Temperature | 5 °C to 55 °C | |
| Protection index | IP 42 | |
| Mounting | Flush mounting (standard box) | |
| Safety | EN 60950-1:2006 | |
| EMC | EN 61000-3-3 | |